



NEONATAL CATHETERIZATION


SPECIFIC INTRODUCERS AND CATHETERS



						PUNCTURE NEEDLE	PEELABLE SHEATH
GAUGE	19	24	24	20	20	24	20
Ø ext (mm)	1 mm	0.7 mm	0.7 mm	0.95 mm	1.1 mm	0.55 mm	1.1 mm
Ø int (mm)	0.86 mm	0.5 mm	0.4 mm	0.66 mm	0.7 mm	0.4 mm	0.7 mm
Compatible catheter size	2 Fr with easylock	1 Fr	1 Fr	2 Fr	1 Fr & 2 Fr	1 Fr & 2 Fr with MST	
Compatible with	epicutan eo cava	premi cath	premi cath	nutriline nutriline twinflo epicutan eo cava	premi cath nutriline nutriline twinflo epicutan eo cava	premi cath nutriline nutriline twinflo epicutan eo cava	

SILICONE CATHETERS


epicutan**eo** cava 2 Fr



1 - 4 kg 2 mm
Ø 0,6mm

POLYURETHANE CATHETERS

premi**cath** 1 Fr



<1 kg 1 mm
Ø 0,35mm

nutriline 2 Fr

Available in 3 and 4 Fr



1 - 4 kg 2 mm
Ø 0,6mm

nutriline twinflo 2 Fr (2 lumens)



1 - 4 kg 2 mm
Ø 0,6mm



Newborn weight



Catheter diameter



Recommended vessel diameter

RECOMMENDATIONS FOR CATHETERIZATION

- 1** Use **sets for neonatal PICC placement**. These sets include specific tools and increase safety. This reduces the risk of infection and ensures the sterility of the process. Sets also save time and are cost-effective.




- 2** **Polyurethane** is more recommended than silicone due to its strength, which facilitates insertion and reduces the risk of complications.


- 3** The size of the catheters should not exceed 1/3 of the internal diameter of the vein. If this is not possible, the catheter should never occupy 45% of the vein's lumen.



CEPHALIC VEIN^{1,2,3,4}
0,8 - 2 mm diameter



BRACHIAL VEIN^{1,2,3,4}
1 - 2,5 mm diameter



BASILIC VEIN^{1,2,3,4}
1,2 - 3 mm diameter

Recommended catheter size

1 Fr

2 Fr

3 Fr

Vein diameter

0.74 mm

1.48 mm

2.22 mm


1 mm

2 mm

3 mm

- 45% - maximum vein/catheter ratio not to be exceeded
- 33% - ideal vein/catheter ratio

- 4** The use of the **Modified Seldinger Technique** (MST) with specific devices such as **microsite®**, presents clear benefits in the PICCs cannulation in neonates. It is a less invasive technique that provides a higher success rate on the first attempt of cannulation.



- 5** Use sutureless fixation systems such as **griplok®**. It is also recommended to use cyanoacrylate glue such as **SecurePortIV®**, at the insertion site. It avoids accidental displacement of the catheter.
- ⚠ **Don't use glue on 2 Fr silicone catheters.**



- 6** Using ultrasound for vein assessment of the patient will allow to choose the most appropriate vein.

- Allow to assess vein size
- Improve the accurate localization of the vessel
- Avoid arterial and nerve punctures
- Increase the success rate on the first attempt
- Reduce the number of punctures
- Decrease patient discomfort
- Minimize the risks of complications associated with the cannulation process

Selecting the correct device ensures safe and efficient venous access, preserving the neonatal venous capital.