WHAT CHARACTERISTICS SHOULD THE IDEAL CONNECTOR FOR NEONATAL ENTERAL NUTRITION HAVE?




SAFE

- Incompatible with Luer connectors and other small diameter connectors used in other applications
- Enables a reliable connection between components
 without leaks or cracks



ADAPTED FOR PREMATURES

- With the smallest possible connectors
- With the lightest possible connectors
- In premature infants, a large and heavy connection carries risks of skin injuries or untimely tube displacement



ACCURATE

- Internal dead space less than 0.1 mL
- Generation of fluid displacement less than 0.05 mL during connection

