
WHAT CHARACTERISTICS SHOULD THE IDEAL CONNECTOR FOR NEONATAL ENTERAL NUTRITION HAVE?



SAFE

 Incompatible with Luer connectors and other small diameter connectors used in other applications

 Enables a reliable connection between components
 without leaks or cracks




ADAPTED FOR PREMATURES

With the smallest possible connectors

With the lightest possible connectors

In premature infants, a
 large and heavy connection
 carries risks of skin injuries
 or untimely tube displacement



ACCURATE

- Internal dead space less than 0.1 mL
- Generation of fluid displacement less than 0.05
 mL during connection

